



Soprano ICE Pretty.Cool.



- The complete laser hair removal solution
- Diode Alex technology - for the widest range of hair types and color
- Speed™ technology - half the treatment time
- Compact technology - versatile treatment methods
- YAG technology - for darker skin types
- Ice™ technology- simultaneous contact cooling



Pixel CO₂ Fractional Ablative CO₂ Laser

- Widest array of handpieces
- World's most advanced scanner
- Fastest treatment



Harmony XL Multi-Application Technology

- Supports Laser, Ultrasound, AFT, NIR and LED technologies
- Over 65 FDA cleared indications
- ClearLift – the first fractional non-ablative Q-switched laser



Legato II Total Scar Solution

- Treat scars, stretch marks and various skin imperfections
- IMPACT trans-epidermal delivery
- Fractional Microplasma technology, ablates and heats the skin through controlled and focused delivery of RF



Accent Ultra V Your New Life Lines

- Face and body contouring, skin tightening and cellulite reduction
- Patented Unipolar technology allows for focused, safe and deep heating
- Cold Ultrasound Shear Wave (™) technology
- Trans-Epidermal Delivery with IMPACT (™)

Alma Lasers is a global innovator of laser, light-based, radiofrequency and ultrasound solutions for the aesthetic and surgical markets. We enable practitioners to offer safe and effective procedures while allowing patients to benefit from state-of-the-art, clinically proven technologies and treatments.

* Design subject to change

Alma Lasers GmbH
Nordostpark 100-102
90411 Nuremberg, Germany
Tel. + 49 911 / 89 11 29-0
Fax + 49 911 / 89 11 29-99
Email: info@almalasers.com
www.almalasers.com
PBAL23121301_01

© 2016 Alma Lasers, Ltd. All rights reserved. Alma Lasers Ltd., its logos are trademarks or registered trademarks of Alma Lasers, Ltd. Product specifications are subject to change without notice.

Designed by **Three Bears Creative Agency**



CONNECT WITH
ALMA LASERS



Alma
Lasers™

True Innovation
the Alma way

THE LARGEST PORTFOLIO



Sopranoice PLATINUM



- The complete laser hair removal
- Trio clustered diode technology
- Combines 3 laser wavelengths into a single innovative handpiece: 755, 810, 1064
- Simultaneously targeting different tissue depths



Alma Q Maximum Power Triple Mode

- Q-switched, Long Pulse Nd: YAG laser, Quasi-LP
- Consecutive pulse technology delivering maximum power
- Pixel fractional Q-Switched Laser
- Unique depth control capabilities
- Safe & effective for all skin types



SINON Pigmented Lesions and Multi-Color Tattoo Removal

- The Q-switched Ruby Laser system with an optimal 694nm wavelength
- Fast treatment - high pulse repetition rate
- Suitable for all skin & types
- New fractional handpiece



Harmony XL Pro Multi-Application Platform

- Combined technologies for better clinical results
- Over 65 FDA - cleared indications
- ClearLift 4D - control the precise depth of the treatment
- ClearSkin - cooled vacuum-assisted laser to treat acne vulgaris



ARION High-end Solution for Hair Removal

- Alexandrite laser system for effective hair removal
- Treatment of vascular lesions and pigmented skin changes
- Burst mode with short sub-pulses for maximum protection of the skin
- Scanner for large-area treatments



308 EXCIMER SYSTEM Small Excimer Impressive Effects


- Monochromatic UVB light therapy
- Clinically proven effectiveness for psoriasis and vitiligo
- Small, mobile and easy to use
- Cost-effective excimer technology



Accent Prime Love Your Shape



- The most advanced platform for body contouring & skin tightening
- Guided ultrasonic wave technology
- Significantly reduced treatment time
- New UltraSpeed Applicator - an extra large applicator plate offers a large coverage rate for optimal results



IDAS Green Light for Mobile Laser Power

- Lbo laser system with 532 nm wavelength
- Vascular indications and pigmented skin changes
- 8W power for effective treatments
- Robust technology with transportation allowance



BURANE II / BURANE II XL Triple Mode - the Simple Choice

- Er: YAG laser system for classical ablative applications
- Burane II XL with triple mode technology: ablation mode, thermal mode and combination mode
- Skin resurfacing, skin rejuvenation and smoothing of scars
- Burane II XL can be used for fractional treatments